


# Nephroline™

Vascular Access Devices

by Meditechlab



**Nephroline™** dialysis catheters cover a wide range of patients' needs.

These catheters are manufactured using the ultimate quality of raw material insuring the physician prescription and the patient's comfort.

The highest technology is mastered to provide:

- Biocompatible polyurethane lumens representing low hypersensitivity risk
- High blood flow due to excellent extrusion
- Kink resistance
- Rotating wing hub for easy suture
- Radio-opaque materials to enhance positioning ease
- Color coded clamps
- Latex-free products

**Nephroline™** Double Lumen catheter kits come on a tray with a catheter of the required specifications (see below), an introducer needle, a wire guide, a dilator and 2 latex-free injection sites/plugs.

## Single Lumen Catheters kit \*

Reference	French size	φ=mm	Length (mm)
SL 8100	8,0	2,7	100
SL 8125	8,0	2,7	125
SL 8150	8,0	2,7	150
SL 8175	8,0	2,7	175
SL 8200	8,0	2,7	200
SL 8225	8,0	2,7	225
SL 8250	8,0	2,7	250

## Double Lumen Catheters kit \*

Reference	Type	French size	φ=mm	Length (mm)
DL 6575	S & J	6,5	2,0	75
DL 65100	S & J	6,5	2,0	100
DL 875	S & J	8,0	2,7	75
DL 8100	S & J	8,0	2,7	100
DL 11100	S & J	11	3,7	100
DL 11125	S & J	11	3,7	125
DL 11150	S & J	11	3,7	150
DL 11175	S & J	11	3,7	175
DL 11200	S & J	11	3,7	200

## Triple Lumen Catheters kit \*

Reference	Type	French size	φ=mm	Length (mm)
TL 12125	S & J	12	4,0	125
TL 12150	S & J	12	4,0	150
TL 12175	S & J	12	4,0	175
TL 12200	S & J	12	4,0	200
TL 12225	S & J	12	4,0	225
TL 12250	S & J	12	4,0	250

\*use **S** for straight (sub-clavian) and **J** for curved (Jugular)

Individual sterile components are available.

Additional references and accessories (such as blades) are available upon request.

